

DISINFECTANT WIPE

Our best-selling range of everyday disinfecting wipes, Medipal Disinfectant Wipes offer great flexibility. Powerful combination chemistry delivers broad spectrum efficacy that's intended for surfaces of non-invasive medical devices.





AT A GLANCE:

- ✓ Alcohol Free
- ✓ PHMB Free
- ✓ Dermatologically Tested

INTENDED FOR:

- ✓ Disinfecting patient shared equipment
- ✓ Use where there is a high risk of infection
- ✓ Disinfecting non-invasive medical devices prior to sterilisation





Scan for online factsheet
and safety data sheet

HOW TO USE:

STEP 1 RISK ASSESSMENT

Please follow your agreed risk assessment policy guidelines regarding the use of PPE. Check device manufacturers for quarternary ammonium compound suitability.



STEP 2 WIPE SELECTION

Check pack for use period and dispense a wipe. Close pack to protect remaining wipes. If visible soil is present, remove with a wipe. Use fresh wipe to disinfect the surface.

STEP 3 S-SHAPE TECHNIQUE


Wipe the surface in an S Shape, moving from clean to dirty. Use the wipe flat not scrunched, ensuring the entire surface is wetted.

Do not go over the same area twice with the same wipe. Use a fresh wipe if your wipe becomes soiled or dry.



STEP 4 DISCARD

Do not re-use wipes. Discard used wipes in the appropriate waste bin following your local agreed guidelines.



Do not flush



STEP 5 LET DRY

Allow the surface to dry naturally before use. Leave the surface for the stated contact time.

Use biocides safely, Always read the label and product information before use. Re-use wipes increases the risk of infection, discard each wipe after single use. Store in cool, dry place out of direct sunlight.

TECHNICAL DETAILS



Medipal's proprietary Tritex[®] material is made by bonding layers of spunbond and meltblown polypropylene fibres together in a high heat process. The addition of a hydrophilic coating means optimum and consistent liquid discharge.

CHEMICAL FORMULATION

Medipal Disinfectant Wipes use a high specification combination chemistry, effectively blending two disinfectant compounds, DDAC (a fourth generation Quat) and Triamine, to produce a highly efficacious, broad spectrum, alcohol free disinfectant.

MATERIALS

Tritex[®] material is made by bonding layers of spunbond and meltblown polypropylene fibres together in a high heat process before a hydrophilic coating is applied. This material is 100% synthetic.



High Strength
The material offers high strength in both directions – across and along the wipe.



Excellent Solution Retention
Ensuring even wetness throughout the wipe.



Very Low Linting
Low risk of leaving contaminating fibres on surfaces.



Efficient Solution Release
Active ingredients are transferred to the surface rather than being trapped in the wipe material fibres.

Efficacy details

Effective Against	Test	Kill Time
BACTERICIDAL		
Acinetobacter baumannii NCTC 12156	EN 13727	30 seconds
E. coli (EBSL) NCTC 10418	EN 13727	30 seconds
E. faecium NCTC 12201	EN 13727	30 seconds
E. hirae NCTC 13383	EN 13727	30 seconds
Enterococcus faecalis (VRE)	EN 13727	1 minute
	EN 13697	1 minute
Enterococcus hirae	EN 16615	1 minute
	EN 13697	1 minute
	EN 14561	1 minute
Escherichia coli	EN 13697	1 minute
K. pneumoniae (CPE) NCTC 7427	EN 13727	30 seconds
P. aeruginosa NCTC 13359	EN 13727	30 seconds
Pseudomonas aeruginosa	EN 16615	1 minute
	EN 13697	1 minute
	EN 14561	1 minute
S. aureus (MRSA)	EN 13727	1 minute
	EN 13697	1 minute
S. aureus NCTC 107788	EN 13727	30 seconds
Staphylococcus aureus	EN 16615	1 minute
	EN 13697	1 minute
	EN 14561	1 minute

Efficacy details

Effective Against	Test	Kill Time
MYCOBACTERICIDAL		
Mycobacterium avium	EN 14348	1 minute
	EN 14563	1 minute
Mycobacterium terrae	EN 14348	1 minute
	EN 14563	1 minute
VIRUCIDAL		
Adenovirus 5*	EN 14476	1 minute
Murine norovirus*	EN 14476	1 minute
Poliovirus*	EN 14476	2 minutes
	*These organisms represent Full Virucidal efficacy according to EN 14476:2013+A2:2019 Annex A	
Vacciniavirus ankara	EN 14476	1 minute
Coronavirus	EN 14476	1 minute
YEASTICIDAL AND FUNGICIDAL		
Candida albicans	EN 16615	1 minute
	EN 13697	1 minute
	EN 13624	1 minute
	EN 14562	2 minutes
Candida auris	EN 13624	1 minute



PRODUCTS WITHIN THIS RANGE

50 WIPE FLOW WRAP



Quantity: 50 per pack
24 packs per case

Classification: Medical Device Class IIa

Product codes include: S524110MPCE S520110MPCE

200 WIPE FLOW WRAP



Quantity: 200 per pack
6 packs per case

Classification: Medical Device Class IIa

Product codes include: S522110MPCE S529110MPCE S625110MPCE

CANISTER



Quantity: 200/240 per pack
10 packs per case

Classification: Medical Device Class IIa

Product codes include: W615110MPCE W520110MPCE

BUCKET



Quantity: 200/225/350 per pack
4 packs per case

Classification: Medical Device Class IIa

Product codes include: W621110MPCE W527110MPCE W522110MPCE

REFILL



Quantity: 225 per pack
4 packs per case

Classification: Medical Device Class IIa

Product codes include: Q527110MPCE

Medipal Disinfectant Wipes Material Compatibility

- ✓ Acrylic (PMMA)
- ✓ Polycarbonate (PC)
- ✓ Polyethylene (PE)
- ✓ Polyvinylchloride (PVC)
- ✓ Waterproof fabric
- ✓ 316 Stainless Steel
- ✓ 304S15 Stainless Steel

