

IPA

Sterile Isopropyl Alcohol 70% v/v with WFI Quality Water

Rapid Action, Low Residue, Fast Drying, Alcohol Disinfectant



USE INSTRUCTIONS:

Apply to surface to ensure complete coverage. Hold spray approximately 15cm to 20cm from area to be treated.

TECHNICAL DATA:

Appearance: Clear liquid free from particles. Specific Gravity (20°C): 0.872 to 0.883. Alcohol proportion: 70% +/- 2% v/v.

STORAGE:

Store upright in original closed containers away from sunlight and extremes of temperature. Full guidance on the handling and disposal of this product is available in the Material Safety Data Sheet (MSDS).

SAFE HANDLING:

Always wear gloves and goggles or face protection. Always read the label and MSDS before use.

MATERIAL COMPATIBILITY:

Application of solutions, when used as directed, will not affect materials normally encountered in the cleanroom. See compatibility information in the Technical File.

BENEFITS:

- Low residue
- Rapid action
- 5L BAG formats manufactured with WFI quality water
- Manufactured to GMP
- Bactericidal and Fungicidal
- Ideal surface and aseptic transfer disinfectant
- Trigger with 1 meter tube-length

TECHNICAL USER INFORMATION:

InSpec IPA 5L BAG is a sterile ready to use cleaner. This is a 70% v/v alcohol solution in WFI quality water.

InSpec IPA is manufactured to GMP in an ISO 5 cleanroom. The product is filled and capped in an ISO class 5 environment.

The solution is filtered through a 0.2 micron filter and gamma irradiated to a validated process (at 25-45kGy) to give a sterility assurance level SAL of 10⁻⁶.

InSpec IPA is supplied with Certificates of Analysis.

All formats of IPA have a stabilised 24-month shelf life.

Microbiological Minimum Efficacy Testing Contact

EN 1276 Bacteria	5 minutes	EN 13697 Bacteria	5 minutes
EN 1650 Fungi	15 minutes	EN 13697 Fungi	15 minutes

InSpec™ IPA Product

Product	Pack Size	Case Size	Product Code
InSpec™ IPA 5L BAG	5L BAG	2 x 5L	IPAWFIBAG-5LS

