

## Tork Blue Couch Roll Advanced



Article	152250
Roll Length	54.5 m
System	C1 - Couch roll system
Roll Width	50 cm
Roll Diameter	12 cm
Number of Sheets	165
Sheet Length	33 cm
Core Inner Diameter	5.1 cm
Ply	2
Embossing	Yes
Colour	Blue

Protect your patient couch with the disposable Tork Couch Roll Advanced. The perforated couch roll is ideal to maintain a high level of hygiene in the healthcare environment. Rolls can be used with the Tork Couch Roll Dispenser.

[www.tork.co.uk](http://www.tork.co.uk)

Better resistance and higher comfort due to the double ply  
Perforation for easy tear-off and consumption control

## Shipping data

	Consumer unit	Transport unit	Pallet
<b>EAN</b>	7310791230375	7310791230320	7322540716146
<b>Pieces</b>	1	9	162
<b>Consumer units</b>	-	9	162
<b>Height</b>	500 mm	500 mm	1161 mm
<b>Width</b>	120 mm	360 mm	1000 mm
<b>Length</b>	120 mm	360 mm	1200 mm
<b>Volume</b>	7.2 dm <sup>3</sup>	64.8 dm <sup>3</sup>	1.2 m <sup>3</sup>
<b>Net weight</b>	980 g	8.8 kg	158.8 kg
<b>Gross weight</b>	1025 g	9.3 kg	167.5 kg
<b>Packing material</b>	Plastic	Plastic	-

## Choose a dispenser



**129184**

**129183**

## Alternative products



**150250**

**150226**

**152226**

## Product certificates



The mark of responsible forestry

## Contact

Essity UK Ltd  
Southfields Road  
Dunstable  
Bedfordshire  
LU6 3EJ

[TorkCS.UK@essity.com](mailto:TorkCS.UK@essity.com)

Tel 01582 677570

Fax 01582 677580