

Tork Couch Roll Advanced



Article	150250
Roll Length	56.1 m
System	C1 - Couch roll system
Roll Width	48 cm
Roll Diameter	12 cm
Number of Sheets	165
Sheet Length	34 cm
Core Inner Diameter	4.4 cm
Ply	2
Embossing	Yes
Colour	White

Protect your patient couch with the disposable Tork Couch Roll Advanced. The perforated couch roll is ideal to maintain a high level of hygiene in the healthcare environment. Rolls can be used with the Tork Couch Roll Dispenser.

www.tork.co.uk

Better resistance and higher comfort due to the double ply

Perforation for easy tear-off and consumption control

Individual plastic wrap for each roll - extra hygienic

Shipping data

	Consumer unit	Transport unit	Pallet
EAN	7310791230351	7310791230344	7322540904192
Pieces	1	9	162
Consumer units	-	9	162
Height	480 mm	480 mm	1121 mm
Width	120 mm	360 mm	1000 mm
Length	120 mm	360 mm	1200 mm
Volume	6.9 dm ³	62.2 dm ³	1.1 m ³
Net weight	969 g	8.7 kg	157 kg
Gross weight	1005 g	9.1 kg	164.3 kg
Packing material	Plastic	Plastic	-

Choose a dispenser



129183



129184

Alternative products



152250



150225



152225

Product certificates



FSC® C030255
The mark of responsible forestry

Contact

Essity UK Ltd

www.tork.co.uk